


# PosiSep® C/C<sup>2</sup> NASAL DRESSINGS

## Engineered for Clinic/Office based nasal procedures

PosiSep® C/C<sup>2</sup> are specifically engineered for use in clinic and office based procedures and are intended to minimize bleeding and edema, prevent adhesions between the septum and the nasal cavity, and act as an adjunct to aid in the natural healing process. PosiSep® C/C<sup>2</sup> offer a slender, easily placed, nasal dressing that allow for gentle placement and rapid expansion upon hydration for minimal bleeding procedures. Both are constructed of a patient-comfortable sponge which is manufactured from naturally occurring chitosan polymers. PosiSep® C<sup>2</sup> expands to a greater volume for patients with no middle turbinate and/or a large middle meatus.

### Designed with patient comfort in mind:

- Gentle streamline structure
- Easy to place
- Utilized in clinic/office
- Great visualization
- Easily conforms to the nasal cavity
- No special storage required
- Biocompatible - made of naturally occurring chitosan polymers
- Rapid expansion upon hydration



### ORDERING INFORMATION

Part Number:	Description:	How Supplied:
9210754	PosiSep® C Nasal Dressing	PosiSep® C & C <sup>2</sup> are packaged five (5) sterile, single use units per carton. One (1) sterile unit contains one (1) device. PosiSep® C & C <sup>2</sup> are terminally sterilized by gamma irradiation.
9210764	PosiSep® C <sup>2</sup> Nasal Dressing	

CAUTION: Federal (U.S.) law restricts this device to sale by or on the order of a physician. PosiSep® C/C<sup>2</sup> is contraindicated for patients who have know allergies to shellfish. For detailed information regarding indications for use, warnings and precautions, see Instructions for Use

To place an order or for more information on PosiSep® C/C<sup>2</sup> or our other medical technologies, please contact a sales representative by calling 866-612-2568.