



# flexibility

at your fingertips.

**NexStat® Topical Hemostat Powder** is intended for use under the care of a health care professional as a topical dressing for the temporary treatment of moderate to severely bleeding wounds such as surgical wounds (post-operative, donor sites, dermatological), cuts and lacerations and for the treatment of mild bleeding from topical ENT surgical wounds and nosebleeds. It is also indicated for control of bleeding from the skin at percutaneous needle access, vascular access and percutaneous catheter access sites.

- Easy to use applicator allows precise application for hard-to-reach surgical wounds
- Biocompatible - made from 100% plant-based polysaccharides
- Versatile - multiple flex tips provided for each application
- No mixing or additional prep time necessary
- No special storage conditions required

**Rx Only.** For detailed information regarding indications for use, warnings and precautions, see Instructions for Use.

### ORDERING INFORMATION

**Part Number:** 9130144      **Description:** NexStat® Topical Hemostat Powder

NexStat® is packaged five (5) sterile units per carton. One (1) sterile unit contains 3 grams of hemostat powder in a compressible bellows and three (3) flexible applicator tips. NexStat® is terminally sterilized by gamma irradiation.



**To place an order or for more information on NexStat® or our other medical technologies, please contact a sales representative by calling 866-612-2568.**