


# NEXFOAM<sup>®</sup> ENT TOPICAL HEMOSTAT SPONGE

bleeding control that's  
**versatile**  
and **powerful.**



**NexFoam<sup>®</sup> Topical Hemostat Sponge** is intended for use under the care of a health care professional as a topical dressing for the temporary treatment of moderate to severely bleeding wounds such as surgical wounds (post-operative, donor sites, dermatological), cuts and lacerations and for the treatment of mild bleeding from topical ENT surgical wounds and nosebleeds. It is also indicated for control of bleeding from the skin at percutaneous needle access, vascular access and percutaneous catheter access sites.

- Can be trimmed or shaped for a variety of applications
- Biocompatible – made from 100% plant-based polysaccharides
- Dissolves naturally over 7-14 days or can be easily removed with gentle irrigation and aspiration
- Easy to apply – no mixing or additional prep time necessary
- No special storage conditions required



**Rx Only.** For detailed information regarding indications for use, warnings and precautions, see Instructions for Use.



## ORDERING INFORMATION

**Part Number:** 9110254

**Description:** NexFoam<sup>®</sup> Topical Hemostatic Sponge (approx. size: 0.6" x 1.9")

NexFoam<sup>®</sup> is packaged five (5) sterile units per carton. One (1) sterile unit contains one (1) hemostatic sponge.  
NexFoam<sup>®</sup> is terminally sterilized by gamma irradiation.

**To place an order or for more information on NexFoam<sup>®</sup> or our other medical technologies, please contact a sales representative by calling 866-612-2568.**