

NEXSTAT[®] PLUS

TOPICAL HEMOSTAT POWDER



flexible control

specially formulated for
ultra absorbency

NexStat[®] Plus Topical Hemostat Powder is intended for use under the care of a health care professional as a topical dressing for the temporary treatment of moderate to severely bleeding wounds such as surgical wounds (post-operative, donor sites, dermatological), cuts and lacerations and for the treatment of mild bleeding from topical ENT surgical wounds and nosebleeds. It is also indicated for control of bleeding from the skin at percutaneous needle access, vascular access and percutaneous catheter access sites.

- Easy to use applicator allows precise application for hard-to-reach surgical wounds
- Biocompatible – made from 100% plant-based polysaccharides
- Versatile – multiple flex tips provided for each application
- No mixing or additional prep time necessary
- No special storage conditions required

New & Improved Formula

NexStat Plus Topical Hemostat Powder is specially formulated with more hydrophillic particles that absorb up to five times their volume. The ultra absorbent particles quickly dehydrate blood cells, resulting in hemoconcentration of platelets and serum proteins to limit and control bleeding.

ORDERING INFORMATION

Part Number:	Description:
9430014	NexStat Plus Topical Hemostat Powder with Applicator (2 grams)
9430034	NexStat Plus Topical Hemostat Powder with Applicator (3 grams)

NexStat Plus Topical Hemostat Powder is packaged with five (5) sterile units per carton. One (1) sterile unit contains 2 grams or 3 grams of hemostat powder in a compressible bellows and three (3) flexible applicator tips. NexStat Plus Topical Hemostat Powder is terminally sterilized by gamma irradiation.

CAUTION: Federal (U.S.) law restricts this device to sale by or on the order of a physician. NexStat is contraindicated for patients who have known allergies to potato starch. For detailed information regarding indications for use, warnings and precautions, see Instructions for Use.



To place an order or for more information on NexStat Plus or our other medical technologies, please contact a sales representative by calling 866-612-2568.